

Anti-Armc6 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Nxe6

Item Number: CM0017902

Introduction

Rabbit polyclonal antibody targeting human ARMC6 (Armadillo repeat-containing protein 6), recommended for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. 54 kDa protein. Ideal for proteomics and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	ARMC6 Rabbit pAb
Remarks/Alias	R30923_1; ARMC6
Species	Human
GeneID (Human)	93436
GeneID	93436
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ARMC6 (NP_001186125.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6NXE6
Protein Weight	54kDa
Shipping	Ice bag