

Anti-Armet/Arp/Manf Rabbit Monoclonal Antibody For Wb, Elisa - P55145

Item Number: CM0002858



Introduction

High-quality rabbit monoclonal antibody against human ARMET/ARP/MANF, validated for WB and ELISA. For neuroscience and ER stress research. Cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	ARMET/ARP/MANF Rabbit mAb
Remarks/Alias	ARP; ARMET; ARMET/ARP/MANF
Species	Human
GeneID (Human)	7873
GeneID	7873
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-182 of human MANF (NP_006001.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55145
Protein Weight	21 kDa
Shipping	Ice bag