

Anti-Armet/Arp/Manf Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P55145

Item Number: CM0025156



Introduction

Rabbit polyclonal antibody against human MANF (ARMET/ARP), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat. Suitable for neurotrophic factor and ER stress studies.

[Learn More](#)

Parameter	Specification
Product Name	ARMET/ARP/MANF Rabbit pAb
Remarks/Alias	ARP; ARMET; ARMET/ARP/MANF
Species	Human
GeneID (Human)	7873
GeneID	7873
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-182 of human ARMET/ARP/MANF (NP_006001.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55145
Protein Weight	21kDa
Shipping	Ice bag