

Anti-Arnt2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hbz2

Item Number: CM0029319

Introduction

Rabbit polyclonal antibody targeting human ARNT2 (bHLHe1) for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying transcription regulation, hypothalamic-pituitary axis development, and nuclear signaling. ~79 kDa.

[Learn More](#)

Parameter	Value
Product Name	ARNT2 Rabbit pAb
Remarks/Alias	WEDAS; bHLHe1; ARNT2
Species	Human
Gene ID	9915
Immunogen	Recombinant fusion protein containing amino acids 440-690 of human ARNT2 (NP_055677.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9HBZ2
Protein Weight	79kDa
Shipping	Ice bag