


Anti-Aromatase (Cyp19A1) Polyclonal Antibody For Wb - P11511

Item Number: CM0024066



Introduction

Rabbit polyclonal antibody targeting human Aromatase (CYP19A1), suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying estrogen biosynthesis and hormone-dependent pathways.

[Learn More](#)

Parameter	Value
Product Name	Aromatase (CYP19A1) Rabbit pAb
Remarks/Alias	ARO; ARO1; CPV1; CYAR; CYP19; CYPXIX; P-450AROM; Aromatase (CYP19A1)
Species	Human
Gene ID (Human)	1588
Gene ID	1588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-503 of human Aromatase (CYP19A1) (NP_000094.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11511
Protein Weight	58kDa
Shipping	Ice bag