


Anti-Aromatase (Cyp19A1) Rabbit Monoclonal Antibody For Wb - P11511

Item Number: CM0002470



Introduction

Rabbit monoclonal antibody against Aromatase (CYP19A1), immunogen: synthetic peptide. Validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Essential for estrogen biosynthesis studies and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Aromatase (CYP19A1) Rabbit mAb
Remarks/Alias	ARO; ARO1; CPV1; CYAR; CYP19; CYPXIX; P-450AROM; Aromatase (CYP19A1)
Species	Human
GeneID (Human)	1588
GeneID	1588
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aromatase (CYP19A1) (P11511).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11511
Protein Weight	58kDa
Shipping	Ice bag