


Anti-Arp2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61160

Item Number: CM0030776



Introduction

Rabbit polyclonal antibody against human Arp2, validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Targets the Arp2/3 complex component involved in actin dynamics and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	Arp2 Rabbit pAb
Remarks/Alias	ARP2; Arp2
Species	Human
GeneID (Human)	10097
GeneID	10097
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human Arp2 (NP_005713.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61160
Protein Weight	45kDa
Shipping	Ice bag