


Anti-Arp2 Rabbit Monoclonal Antibody For Wb And Elisa - P61160

Item Number: CM0010043



Introduction

Rabbit monoclonal antibody against human Arp2 (ARP2), validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying actin polymerization, cell motility, and nuclear DNA repair.

[Learn More](#)

Parameter	Value
Product Name	Arp2 Rabbit mAb
Remarks/Alias	ARP2; Arp2
Species	Human
GeneID (Human)	10097
Immunogen	A synthetic peptide within amino acids 295-394 of human Arp2 (P61160)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61160
Protein Weight	45kDa
Shipping	Ice bag