


Anti-Arp3 Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - P61158

Item Number: CM0011906



Introduction

Rabbit monoclonal anti-Arp3 antibody for Western blot, IHC-P, IF/ICC, and ELISA. Detects human, mouse, rat Arp3 (ACTR3) at 47kDa. Suitable for cytoskeleton, cell motility, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Arp3 Rabbit mAb
Remarks/Alias	ARP3; Arp3
Species	Human
GeneID (Human)	10096
GeneID	10096
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 319-418 of human Arp3 (P61158).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61158
Protein Weight	47kDa
Shipping	Ice bag