


Anti-Arpc2 Monoclonal Antibody For Wb - 015144

Item Number: CM0003102



Introduction

High-quality anti-ARPC2 rabbit monoclonal antibody validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat ARPC2, a 34 kDa subunit of the Arp2/3 complex essential for actin polymerization and cell motility. Ideal for cytoskeleton and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	ARPC2 Rabbit mAb
Remarks/Alias	ARC34; PRO2446; p34-Arc; PNAS-139; ARPC2
Species	Human
GeneID (Human)	10109
GeneID	10109
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15144
Protein Weight	34kDa
Shipping	Ice bag