


Anti-Arpc2 Polyclonal Antibody For Wb, Ihc-P, Elisa - O15144

Item Number: CM0014659



Introduction

Anti-ARPC2 rabbit polyclonal antibody for WB, IHC-P, ELISA. Targets human ARPC2 (p34-Arc), component of Arp2/3 complex, involved in actin polymerization. Validated in Human, Mouse, Rat.

[Learn More](#)

Core Parameter	Specification
Product Name	ARPC2 Rabbit pAb
Remarks/Alias	ARC34; PRO2446; p34-Arc; PNAS-139; ARPC2
Species	Human
Gene ID (Human)	10109
Gene ID	10109
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ARPC2 (NP_690601.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15144
Protein Weight	34kDa
Shipping	Ice bag