


## Anti-Arpc3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15145

Item Number: CM0030860



### Introduction

Rabbit polyclonal antibody against human ARPC3, suitable for WB, IF/ICC, and ELISA. Detects endogenous levels of actin-related protein 2/3 complex subunit 3 (p21-ARC) in human, mouse, and rat samples. Ideal for studying actin polymerization and DNA repair.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | ARPC3 Rabbit pAb                                                                                                    |
| Remarks/Alias    | ARC21; p21-Arc; ARPC3                                                                                               |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 10094                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human ARPC3 (NP_001265485.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | O15145                                                                                                              |
| Protein Weight   | 21kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |