


Anti-Arpc5 (P16-Arc) Rabbit Monoclonal Antibody For Wb And Elisa - O15511

Item Number: CM0012716



Introduction

Rabbit monoclonal antibody targeting human ARPC5/p16-ARC, a subunit of the Arp2/3 complex involved in actin polymerization and DNA repair. Validated for WB and ELISA in human, mouse, rat. Suitable for cell motility and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	p16-ARC/ARC16/ARPC5 Rabbit mAb
Remarks/Alias	ARC16; p16-Arc; dj127C7.3; p16-ARC/ARC16/ARPC5
Species	Human
GeneID (Human)	10092
GeneID	10092
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human p16-ARC/ARC16/ARPC5 (O15511)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15511
Protein Weight	16kDa
Shipping	Ice bag