

Anti-Arrb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P32121

Item Number: CM0020729



Introduction

High-quality rabbit polyclonal antibody targeting human ARR2 (beta-arrestin-2), suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying GPCR desensitization, endocytosis, and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	ARRB2 Rabbit pAb
Remarks/Alias	ARB2; ARR2; BARR2; ARRB2
Species	Human
GeneID (Human)	409
GeneID	409
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 295-394 of human ARRB2 (NP_945355.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32121
Protein Weight	46kDa
Shipping	Ice bag