

Anti-Arrdc4 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nct1

Item Number: CM0020267

Introduction

Rabbit polyclonal antibody targeting human ARRDC4, validated for WB and ELISA. Suitable for studying arrestin domain-containing protein 4 functions in GPCR endocytosis, glucose uptake, and innate immune response. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	ARRDC4 Rabbit pAb
Remarks/Alias	ARRDC4
Species	Human
Gene ID	91947
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ARRDC4 (NP_899232.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NCT1
Protein Weight	45kDa
Shipping	Ice bag