


Anti-Arsa Polyclonal Antibody For Wb And Elisa - P15289

Item Number: CM0023390



Introduction

Rabbit polyclonal antibody against human ARSA, validated for Western blot (WB) and ELISA applications. Cross-reacts with mouse. Ideal for research on lysosomal sulfatide metabolism and metachromatic leukodystrophy.

[Learn More](#)

Parameter	Value
Product Name	ARSA Rabbit pAb
Remarks/Alias	ASA; MLD; ARSA
Species	Human
GeneID (Human)	410
GeneID	410
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-509 of human ARSA (NP_000478.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P15289
Protein Weight	54kDa
Shipping	Ice bag