

Anti-Art1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P52961

Item Number: CM0016706

Introduction

Rabbit polyclonal antibody against human ART1 (ARTC1/CD296) for WB, IHC-P, IF/ICC, ELISA. Detects endogenous ART1 at ~36 kDa in human and mouse; ideal for ADP-ribosylation studies.

[Learn More](#)

Parameter	Value
Product Name	ART1 Rabbit pAb
Remarks/Alias	RT6; ART2; ARTC1; CD296; ART1
Species	Human
Gene ID (Human)	417
Gene ID	417
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-295 of human ART1 (NP_004305.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52961
Protein Weight	36kDa
Shipping	Ice bag