

Anti-Art3 Rabbit Pab For Wb, If-P, Elisa - Q13508

Item Number: CM0024039

Introduction

High-quality rabbit polyclonal antibody against human ART3 (ARTC3), ecto-ADP-ribosyltransferase 3. Validated for WB, IF-P, ELISA; cross-reacts with rat. Suitable for testis-specific protein detection and cell membrane studies.

[Learn More](#)

Parameter	Value
Product Name	ART3 Rabbit pAb
Remarks/Alias	ARTC3; ART3
Species	Human
GeneID (Human)	419
GeneID	419
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human ART3 (NP_001123488)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13508
Protein Weight	44kDa
Shipping	Ice bag