

Anti-Art5 Polyclonal Antibody For Wb, Elisa - Q96L15

Item Number: CM0024938

Introduction

Rabbit polyclonal antibody against human ART5 (ARTC5) validated for Western blotting (WB) and ELISA. Detects human, mouse, and rat ART5 (32 kDa). Suitable for ecto-ADP-ribosyltransferase 5 research.

[Learn More](#)

Parameter	Value
Product Name	ART5 Rabbit pAb
Remarks/Alias	ARTC5; ART5
Species	Human
GeneID (Human)	116969
GeneID	116969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-200 of human ART5 (NP_443750.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96L15
Protein Weight	32kDa
Shipping	Ice bag