


Anti-Artn Rabbit Polyclonal Antibody For Wb, Elisa - Q5T4W7

Item Number: CM0030343



Introduction

Rabbit polyclonal antibody against human Artemin (ARTN) for Western blot and ELISA. Cross-reacts with mouse and rat. Useful for investigating peripheral neuron survival, dopaminergic neuron support, and RET/GFR α 3 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	ARTN Rabbit pAb
Remarks/Alias	ART; EVN; NBN; ENOVIN; ARTN
Species	Human
GeneID (Human)	9048
GeneID	9048
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 122-220 of human ARTN (NP_476432.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5T4W7
Protein Weight	23kDa
Shipping	Ice bag