

Anti-Arvcf Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00192

Item Number: CM0023380

Introduction

Rabbit polyclonal antibody targeting human ARVCF, validated for WB, IF/ICC, ELISA. Detects mouse and rat orthologs. Suitable for research on alternative splicing and cell adhesion mechanisms.

[Learn More](#)

Parameter	Value
Product Name	ARVCF Rabbit pAb
Remarks/Alias	ARVCF
Species	Human
GeneID (human)	421
GeneID	421
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ARVCF (NP_001661.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O00192
Protein Weight	105kDa
Shipping	Ice bag