

Anti-Arx Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Qs3

Item Number: CM0020607

Introduction

Rabbit polyclonal antibody against human ARX, a homeobox transcription factor involved in brain development and pancreatic cell differentiation. Suitable for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	ARX Rabbit pAb
Remarks/Alias	ISSX; PRTS; CT121; EIEE1; MRX29; MRX32; MRX33; MRX36; MRX38; MRX43; MRX54; MRX76; MRX87; MRXS1; ARX
Species	Human
Gene ID (Human)	170302
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human ARX (NP_620689.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Qs3
Protein Weight	58kDa
Shipping	Ice bag