

Anti-Asah1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13510

Item Number: CM0015225



Introduction

Rabbit polyclonal antibody targeting human acid ceramidase (ASA1), offered for Western blot, immunofluorescence/immunocytochemistry (IF/ICC), and ELISA. Suitable for studying sphingolipid metabolism and ceramide-related signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	ASA1 Rabbit pAb
Remarks/Alias	AC, PHP, ASAH, PHP32, ACDase, SMAPME, ASAH1
Species	Human
GeneID (Human)	427
GeneID	427
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-389 of human ASA1 (NP_001120977.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13510
Protein Weight	45kDa
Shipping	Ice bag