

Anti-Asap2 Polyclonal Antibody For Wb - O43150

Item Number: CM0022113

Introduction

Rabbit polyclonal antibody against human ASAP2 (DDEF2/PAG3). Validated for WB and ELISA, with mouse cross-reactivity. Ideal for research in vesicular trafficking and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	ASAP2 Rabbit pAb
Remarks/Alias	PAP; PAG3; AMAP2; DDEF2; SHAG1; CENTB3; Pap-alpha; ASAP2
Species	Human
Gene ID (Human)	8853
Gene ID	8853
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human ASAP2 (NP_001128663)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O43150
Protein Weight	112kDa
Shipping	Ice bag