

Anti-Asb13 Polyclonal Antibody For Wb, Elisa - Q8Wxk3

Item Number: CM0015121

Introduction

Rabbit polyclonal antibody against human ASB13, validated for WB and ELISA. Cross-reacts with mouse and rat. ASB13 is a substrate-recognition component of SCF-like E3 ubiquitin ligase, mediating protein degradation. Ideal for ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	ASB13 Rabbit pAb
Remarks/Alias	ASB13
Species	Human
Gene ID (Human)	79754
Gene ID	79754
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-278 of human ASB13 (NP_078977.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WXX3
Protein Weight	30kDa
Shipping	Ice bag