

Anti-Asb17 Polyclonal Antibody For Wb, Elisa - Q8Vhp9

Item Number: CM0021669

Introduction

Rabbit polyclonal antibody against ASB17 (Ankyrin repeat and SOCS box protein 17). Suitable for WB and ELISA. Cross-reacts with human and rat species. Protein weight approx. 34 kDa. For research use only.

[Learn More](#)

Parameter	Value
Product Name	ASB17 Rabbit pAb
Remarks/Alias	4933436O18Rik; ASB17
Species	Human
GeneID (Human)	66772
GeneID	66772
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 78-177 of mouse ASB17 (NP_080034.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8VHP9
Protein Weight	34kDa
Shipping	Ice bag