

Anti-Asb3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y575

Item Number: CM0023040

Introduction

High-quality rabbit polyclonal antibody against human ASB3, validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat ASB3. Ideal for ubiquitin-proteasome pathway research.

[Learn More](#)

Parameter	Value
Product Name	ASB3 Rabbit pAb
Remarks/Alias	ASB-3; ASB3
Species	Human
GeneID (Human)	51130
GeneID	51130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human ASB3 (NP_057199.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y575
Protein Weight	58kDa
Shipping	Ice bag