

Anti-Asb4 Polyclonal Antibody For Ihc-P, Elisa - Q9Y574

Item Number: CM0024663

Introduction

High-quality anti-ASB4 rabbit polyclonal antibody validated for IHC-P and ELISA. Recognizes human and mouse ASB4, a SOCS-box E3 ligase involved in vascular development and ubiquitination. Ideal for research use.

[Learn More](#)

Parameter	Description
Product Name	ASB4 Rabbit pAb
Remarks/Alias	ASB-4; ASB4
Species	Human
GeneID (Human)	51666
GeneID	51666
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ASB4 (NP_665879.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y574
Protein Weight	48kDa
Shipping	Ice bag