

Anti-Ascc2 Polyclonal Antibody For Wb, Elisa - Q9H1I8

Item Number: CM0013786



Introduction

Rabbit polyclonal antibody against human ASCC2, validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for DNA repair and ribosome quality control research.

[Learn More](#)

Parameter	Value
Product Name	ASCC2 Rabbit pAb
Remarks/Alias	p100; ASC1p100; ASCC2
Species	Human
GeneID (Human)	84164
GeneID	84164
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 488-757 of human ASCC2 (NP_115580.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1I8
Protein Weight	86kDa
Shipping	Ice bag