


Anti-Ascc3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N3C0

Item Number: CM0027892



Introduction

Rabbit polyclonal antibody targeting human ASCC3 (asc1p200/HELIC1), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Detects a ~251 kDa protein involved in DNA alkylation repair and ribosome quality control. Ideal for research use. UniProt Q8N3C0.

[Learn More](#)

Parameter	Specification
Product Name	ASCC3 Rabbit pAb
Remarks/Alias	RNAH; HELIC1; ASC1p200; ASCC3
Species	Human
GeneID (Human)	10973
GeneID	10973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human ASCC3 (NP_071374.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N3C0
Protein Weight	251kDa
Shipping	Ice bag