

Anti-Ascl4 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q6Xd76

Item Number: CM0006492

Introduction

High-quality anti-ASCL4 rabbit monoclonal antibody for WB, IF/ICC, and ELISA. Specific for human ASCL4, a transcription factor involved in skin development. Suitable for dermatology and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	ASCL4 Rabbit mAb
Remarks/Alias	ASH-4; HASH4; bHLHa44; ASCL4
Species	Human
GeneID	121549 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human ASCL4 (NP_982260.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q6XD76
Protein Weight	19kDa
Shipping	Ice bag