

Anti-Ascl4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Xd76

Item Number: CM0016351

Introduction

High-quality rabbit polyclonal antibody targeting human ASCL4 (bHLHa44), validated for Western blot, IHC-P, and ELISA in human, mouse, and rat samples. Suitable for skin development and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	ASCL4 Rabbit pAb
Remarks/Alias	ASH-4; HASH4; bHLHa44; ASCL4
Species	Human
GeneID (Human)	121549
GeneID	121549
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human ASCL4 (NP_982260.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6XD76
Protein Weight	19 kDa
Shipping	Ice bag