


Anti-Asct2/Slc1A5 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15758

Item Number: CM0006849



Introduction

Rabbit monoclonal antibody targeting human ASCT2/SLC1A5 (Neutral amino acid transporter B(0)), validated for WB, IHC-P, ELISA. Ideal for research on glutamine metabolism, D-serine transport, and retroviral entry.

[Learn More](#)

Parameter	Value
Product Name	ASCT2/SLC1A5 Rabbit mAb
Remarks/Alias	R16; AAAT; ATBO; M7V1; RDRC; ASCT2; M7VS1; ASCT2/SLC1A5
Species	Human
Gene ID (Human)	6510
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 483-541 of human ASCT2/SLC1A5(NP_005619.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q15758
Protein Weight	57kDa
Shipping	Ice bag