


Anti-Asc/Tms1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9ULz3

Item Number: CM0007471



Introduction

Rabbit monoclonal antibody against human ASC/TMS1 (PYCARD) for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Key adaptor in apoptosis and inflammasome signaling. Suitable for immunology and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	ASC/TMS1 Rabbit mAb
Remarks/Alias	ASC; TMS; TMS1; CARD5; TMS-1; ASC/TMS1
Species	Human
Gene ID (Human)	29108
Gene ID	29108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 4-110 of human ASC/TMS1 (NP_037390.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULZ3
Protein Weight	22kDa
Shipping	Ice bag