


## Anti-Asc/Tms1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9ULz3

Item Number: CM0023109



### Introduction

Rabbit polyclonal antibody targeting human ASC/TMS1 (Q9ULZ3). Suitable for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for apoptosis and inflammasome research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | ASC/TMS1 Rabbit pAb                                                                                                |
| Remarks/Alias    | ASC; TMS; TMS1; CARD5; TMS-1; ASC/TMS1                                                                             |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 29108                                                                                                              |
| GeneID           | 29108                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-91 of human ASC/TMS1 (NP_037390.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q9ULZ3                                                                                                             |
| Protein Weight   | 22kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |