


Anti-Asgr1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P07306

Item Number: CM0027591



Introduction

Rabbit polyclonal antibody targeting human ASGR1, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. ASGR1 is a hepatocyte-specific lectin mediating glycoprotein clearance, platelet removal, and SARS-CoV-2 entry.

[Learn More](#)

Parameter	Value
Product Name	ASGR1 Rabbit pAb
Remarks/Alias	HL-1; ASGPR; ASGPR1; CLEC4H1; ASGR1
Species	Human
GeneID (Human)	432
GeneID	432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human ASGR1 (NP_001662.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07306
Protein Weight	33kDa
Shipping	Ice bag