


## Anti-Asgr2 Polyclonal Antibody For Wb, Elisa - P07307

Item Number: CM0018299



### Introduction

CamelBio's Anti-ASGR2 polyclonal antibody (rabbit) detects human ASGR2, cross-reacts with mouse and rat, and is validated for WB and ELISA. Ideal for studying endocytosis of asialoglycoproteins.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | ASGR2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | HL-2; HBXBP; ASGPR2; ASGP-R2; CLEC4H2; ASGR2                                                                       |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 433                                                                                                                |
| GeneID           | 433                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 131-287 of human ASGR2 (NP_550435.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Mouse, Rat                                                                                                         |
| SWISS            | P07307                                                                                                             |
| Protein Weight   | 35 kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |