


Anti-Asgr2 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P07307

Item Number: CM0009915



Introduction

Anti-ASGR2 rabbit monoclonal antibody validated for WB, IF-P, ELISA. Targets human asialoglycoprotein receptor 2 (HL-2, CLEC4H2), a hepatic lectin mediating endocytosis of desialylated glycoproteins. For research use.

[Learn More](#)

Parameter	Value
Product Name	ASGR2 Rabbit mAb
Remarks/Alias	HL-2; HBXBP; ASGPR2; ASGP-R2; CLEC4H2
Species	Human
GeneID (Human)	433
GeneID	433
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-311 of human ASGR2 (NP_550435.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P07307
Protein Weight	35kDa
Shipping	Ice bag