


Anti-Asgr2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P07307

Item Number: CM0029411



Introduction

Rabbit polyclonal antibody against human ASGR2 (Asialoglycoprotein receptor 2), validated for WB, IF-P, ELISA; reacts with human and mouse. Ideal for research on hepatic glycoprotein endocytosis and liver targeting.

[Learn More](#)

Parameter	Value
Product Name	ASGR2 Rabbit pAb
Remarks/Alias	HL-2; HBXBP; ASGPR2; ASGP-R2; CLEC4H2; ASGR2
Species	Human
GeneID (Human)	433
GeneID	433
Immunogen	Recombinant protein corresponding to amino acids 58-287 of human ASGR2 (NP_550435.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07307
Protein Weight	35kDa
Shipping	Ice bag