

Anti-Ash2L Monoclonal Antibody For Wb, Ihc, If, Elisa - Q9Ubl3

Item Number: CM0010668



Introduction

Recombinant rabbit monoclonal antibody against human ASH2L (ASH2, Bre2). Recognizes human, mouse, and rat. Suitable for WB, IHC-P, IF/ICC, IF-P, and ELISA. UniProt Q9UBL3.

[Learn More](#)

Parameter	Value
Product Name	ASH2L Rabbit mAb
Remarks/Alias	ASH2, Bre2, ASH2L1, ASH2L2, ASH2L
Species	Human
GeneID (Human)	9070
GeneID	9070
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 529-628 of human ASH2L (Q9UBL3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBL3
Protein Weight	69kDa
Shipping	Ice bag