


Anti-Asl Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04424

Item Number: CM0025272



Introduction

High-quality rabbit polyclonal antibody against human ASL (argininosuccinate lyase), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for urea cycle and arginine metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ASL Rabbit pAb
Remarks/Alias	ASAL; ASL
Species	Human
GeneID (Human)	435
GeneID	435
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 289-464 of human ASL (NP_000039.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04424
Protein Weight	52kDa
Shipping	Ice bag