

Anti-Asmt Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P46597

Item Number: CM0027805



Introduction

Anti-ASMT rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets acetylserotonin O-methyltransferase, a key enzyme in melatonin biosynthesis. Ideal for pineal gland and circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	ASMT Rabbit pAb
Remarks/Alias	ASMTY; HIOMT; HIOMTY; ASMT
Species	Human
Gene ID (Human)	438
Gene ID	438
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-298 of human ASMT (NP_001164510.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46597
Protein Weight	38kDa
Shipping	Ice bag