

Anti-Asmtl Rabbit Polyclonal Antibody For Wb, Ip, Elisa - O95671

Item Number: CM0030186



Introduction

Rabbit polyclonal antibody targeting human ASMTL (ASTML/ASMTLX/ASMTLY), validated for WB, IP, ELISA. Reacts with human and mouse. Targets nucleotide metabolism enzyme dTTP/UTP pyrophosphatase/methyltransferase.

[Learn More](#)

Parameter	Value
Product Name	ASMTL Rabbit pAb
Remarks/Alias	ASTML; ASMTLX; ASMTLY; ASMTL
Species	Human
Gene ID (Human)	8623
Gene ID	8623
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ASMTL (NP_004183.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95671
Protein Weight	69kDa
Shipping	Ice bag