

Anti-Asna1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43681

Item Number: CM0029443



Introduction

Rabbit polyclonal antibody targeting human ASNA1 (GET3/TRC40), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying tail-anchored protein delivery and ER pathways.

[Learn More](#)

Parameter	Value
Product Name	ASNA1 Rabbit pAb
Remarks/Alias	ARSA1; ASNA1; CMD2H; TRC40; ARSA-I; ASNA-I
Species	Human
GeneID (Human)	439
GeneID	439
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-348 of human ASNA1 (NP_004308.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43681
Protein Weight	39kDa
Shipping	Ice bag