


Anti-Asns Rabbit Monoclonal Antibody For Wb, Elisa - P08243

Item Number: CM0002823



Introduction

Anti-ASNS rabbit monoclonal antibody validated for WB and ELISA. Cross-reactive with human, mouse, and rat. Targets asparagine synthetase (P08243), a glutamine-dependent ligase essential for asparagine biosynthesis and cell cycle control.

[Learn More](#)

Parameter	Value
Product Name	ASNS Rabbit mAb
Remarks/Alias	TS11; ASNSD; ASNS
Species	Human
Gene ID (Human)	440
Gene ID	440
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human ASNS (P08243).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P08243
Protein Weight	64kDa
Shipping	Ice bag