

Anti-Asns Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08243

Item Number: CM0026251



Introduction

High-quality anti-ASNS rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying asparagine biosynthesis and cell cycle control.

[Learn More](#)

Parameter	Specification
Product Name	ASNS Rabbit pAb
Remarks/Alias	TS11; ASNSD; ASNS
Species	Human
GeneID (Human)	440
GeneID	440
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 212-561 of human ASNS (NP_001664.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08243
Protein Weight	64kDa
Shipping	Ice bag