

Anti-Asns Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08243

Item Number: CM0016742



Introduction

CamelBio offers anti-ASNS rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets asparagine synthetase (P08243), involved in asparagine biosynthesis and cell cycle control.

[Learn More](#)

Parameter	Value
Product Name	ASNS Rabbit pAb
Remarks/Alias	TS11, ASNSD, ASNS
Species	Human
GeneID (human)	440
GeneID	440
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 262-561 of human ASNS (NP_001664.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08243
Protein Weight	64kDa
Shipping	Ice bag