


Anti-Aspa Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P45381

Item Number: CM0026158



Introduction

Rabbit polyclonal antibody against aspartoacylase (ASPA) for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for studying N-acetylaspartate metabolism, white matter maintenance, and Canavan disease.

[Learn More](#)

Parameter	Value
Product Name	ASPA Rabbit pAb
Remarks/Alias	ASP; ACY2; ASPA
Species	Human
GeneID (Human)	443
GeneID	443
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human ASPA (NP_001121557.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45381
Protein Weight	36kDa
Shipping	Ice bag