


Anti-Asph Polyclonal Antibody For If/Icc, Elisa - Q12797

Item Number: CM0030970



Introduction

This rabbit polyclonal antibody recognizes human ASPH, also known as aspartate beta-hydroxylase. It is suitable for immunofluorescence (IF/ICC) and enzyme-linked immunosorbent assay (ELISA), supporting research into calcium signaling and post-translational modifications.

[Learn More](#)

Parameter	Value
Product Name	ASPH Rabbit pAb
Remarks/Alias	AAH; BAH; HAAH; JCTN; FDLAB; junctin; CASQ2BP1; ASPH
Species	Human
GeneID (Human)	444
GeneID	444
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 81-270 of human ASPH (NP_001158227.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q12797
Protein Weight	86kDa
Shipping	Ice bag