


Anti-Aspscr1 (Tug) Rabbit Polyclonal Antibody For Wb, If/Icc, Ip - Q9Bze9

Item Number: CM0031001



Introduction

Rabbit polyclonal antibody targeting ASPSCR1 (TUG) for Western blot, immunofluorescence, immunoprecipitation, and ELISA. Reacts with human, mouse, rat. Ideal for GLUT4 trafficking and insulin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	TUG/ASPSCR1 Rabbit pAb
Remarks/Alias	TUG; ASPL; ASPS; RCC17; UBXD9; UBXN9; ASPCR1
Species	Human
GeneID (human)	79058
GeneID	79058
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 284-553 of human TUG/ASPSCR1 (NP_076988.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZE9
Protein Weight	60kDa
Shipping	Ice bag